

DISC LIST

Drug Name	Code	Strengths for Standardised Methodologies						OTHER STRENGTHS					
		BSAC			EUCAST		CLSI						
Amikacin	AMK	30			30		30	10	15				
Amoxicillin	AMX	2						5	10	25	30		
Amoxicillin/Clavulanic acid	AUG	3	30		30		30						
Ampicillin	AMP	2	10	25	2	10	10	30					
Ampicillin / sulbactam	UNZ				10/10		10/10						
Ampiclox	ACX							30					
Apramycin	APR							15	30				
Azithromycin	AZM	15					15	125					
Azlocillin	AZL						75	30					
Aztreonam	AZT	30			30		30						
Bacitracin	BAC							10					
Carbenicillin	CAR	100					100						
Cefacetrile	CCT							30					
Cefaclor	CCL	30			30		30						
Cefadroxil	CDR				30								
Cefamandole	CMD	30					30						
Cefatrizine	CTZ							30					
Cefazolin	CZL				30		30						
Cefdinir	CNR						5						
Cefepime	CFP	30			30		30						
Cefepime/Clavulanate	CZP	30/10											
Cefixime	CXM	5			5		5	30					
Cefmetazole	CMZ						30						
Cefodizime	CDZ							30					
Cefoperazone	CPZ	30					75						
Cefoperazone/Sulbactam	CPS							50					
Cefotaxime	CTX	5	30		5		30	25					
Cefotaxime/Clavulanate	CXV						30/10						
Cefotetan	CTT	30					30						
Cefotiam	CTM							30					
Cefoxitin	CXT	10	30		30		30						
Cefpirome	CPM	30											
Cefpodoxime	CPD	1	10		10		10	30					
Cefpodoxime/Clavulanate	CDV	10/1	30/10										
Cefprozil	CFO						30						
Cefsulodin	CDN							30					
Ceftazidime	CAZ	30			10		30						
Ceftazidime/Clavulanate	CZV	30/10					30/10						
Ceftibutin	CTB				30		30	10					
Ceftiofur	CFU							20					
Ceftizoxime	CZM	30											
Ceftrazidime	CTZ							30					
Ceftriaxone	CTR	5	30		30		30						
Cefuroxime	CRX	5	30		30		30						
Cephalexin	CLX	30			30			5					
Cephalothin	CLT						30						
Cephradine	CRD	30						15					
Chloramphenicol	CHL	10	30		30		30	25					
Chlortetracycline	CHT							25	30				
Cinoxin	CNX							100					
Ciprofloxacin	CPR	1	5		5			2	10	30			
Clarithromycin	CLM	2	5				15						
Clindamycin	CLN	2			2		2	10	30				
Cloxacillin	CXC							5	10				
Colistin sulphate	COL	25					10	50	100				
Cotrimoxazole	COT	25			25		25	10	30	100			

Drug Name	Code	Strengths for Standardised Methodologies						OTHER STRENGTHS			
		BSAC		EUCAST		CLSI					
Difloxacin	DIF							25			
Doripenem	DOR			10		10					
Doxycycline	DOX	30						10			
Enrofloxacin	ENO							5	10		
Ertapenem	ERT	10		10		10					
Erythromycin	ERY	5		15		15		2	10		
Flavomycin	FLV							10			
Florfenicol	FFL							30			
Flucloxacillin	FLX							5			
Flumequine	FLQ							4	30		
Fosfamycin	FOS	200				200					
Framycetin	FRA							50	100		
Furazolidone	FUR							10	20	50	100
Fusidic acid	FUS	10		10							
Gatifloxacin	GAT	2				5					
Gemifloxacin	GEM					5					
Gentamicin	GEN	10	200	10	30	10	120	250			
Imipenem	IMP	10		10		10		30			
Isepamicin	ISE							30			
Isoconazole	ISO							50			
Josamycin	JOS							15			
Kanamycin	KAN					30		120	1000		
Levofloxacin	LEV	1	5	5		5					
Lincomycin	LIN							2	10	15	
Linezolid	LNZ	10		10		30					
Lomefloxacin	LOM					10					
Loracarbef	LOR					30		10			
Marbofloxacin	MAR							5			
Mecillinam	MEC	10	50	10		10		25			
Meropenem	MER	10		10		10					
Metronidazole	MTZ	5						50			
Mezlocillin	MEZ	75				75		30			
Minocycline	MIN	30		30		30					
Moxalactam	MOX	1	5			30					
Moxifloxacin	MXF			5		5		1			
Mupirocin	MUP	5	20	200		200					
Nafcillin	NAF					1					
Nalidixic acid	NAL	30				30					
Neomycin	NEO	10						30			
Netilmicin	NET	10		10		30					
Nitrofurantoin	NIT	200		100		300		30	50		
Nitroxilin	NIC				10			30	50		
Norfloxacin	NOR	2		1		10		5			
Novobiocin	NOV							5	10	30	
Ofloxacin	OFL	5		1	5	5		10	25	30	
Oleandomycin	OLN							15	50		
Oxacillin	OXC	1		1		1		10			
Oxolinic acid	OXA							2	4		
Oxytetracycline	OXT							20	25	30	
Pefloxacin	PEF							5	30		

Drug Name	Code	Strengths for Standardised Methodologies			OTHER STRENGTHS					
		BSAC		EUCAST	CLSI					
Penicillin G	PEN	1		1		10	0.03	0.25	2	5
Penicillin-phenoxymethyl	PNV			10						
Pipemidic Acid	PPA						20			
Piperacillin	PIP	75		30		100				
Piperacillin / tazobactam	TZP	75/10		30/8		100/10				
Polymyxin B	POL					300	30	100		
Pristinamycin	PRS						15			
Quinupristin/Dalfopristin	SYN	15		15		15				
Ramoplanin										
Razupenem	RZM					10				
Rifabutin	RFB									
Rifaximin										
Rifampicin	RIF	2	5	5		5	10	15	30	
Rifamycin	RMC						30			
Roxithromycin	ROX						30			
Sparfloxacin	SAR					5				
Sisomicin	SIS						10			
Spectinomycin	SPC	25				100	10	500		
Spiramycin	SPR						30	100		
Streptomycin	STR	10		300		10	30	100		
Sulbactam	SBC						10			
Sulfafurazole	SFZ					250	300	1000		
Sulphamethoxazole	SMX	100								
Sulphonamides	SUL					250	300	200		
Teicoplanin	TCP	30		30		30	5			
Telithromycin	TLM	15		15		15				
Televancin						30				
Temocillin	TEM	30								
Tetracycline	TET	10		30		30				
Thiamphenicol	THI						30			
Tiamulin	TAM						30			
Ticarcillin	TIC	75		75		75	15			
Ticarcillin/clavulanic	BBL	75/10		75/10		75/10				
Tigecycline	TIG	15		15		15				
Tilmicosin	TIL						15	30		
Tobramycin	TOB	10		10		10	30			
Trimethoprim	TRM	2.5	5	5		5	1.25	2.5	25	100
Trimethoprim/Sulphamethoxazole	COT	25		25		25	10	30	100	
Trovafoxacin	TRV					10				
Tylosin	TYL						30			
Unazid	UNZ			10/10		10/10				
Uticillin	UTI						100			
Valnemulin	VAL						30			
Vancomycin	VAN	5		5		30				

CLSI Clinical and Laboratory Standards Institute

www.clsi.org

BSAC British Society for Antimicrobial Chemotherapy

www.bsac.org.uk

EUCAST European Committee on Antimicrobial Susceptibility Testing

www.eucast.org

ANTIFUNGAL DRUGS

Drug Name	Code	Disc Strength		
Amphotericin	AMT	10	20	
Clotrimazole	CLO	10		
Econazole	ECN	10		
Fluconazole	FLC	10	25	
Flucytosine	FLU	1		
Griseofulvin	GFL	10		
Itraconazole	ITR	1	10	
Ketoconazole	KET	1	10	
Miconazole	MIC	5	10	50
Natamycin	NAT	10		
Nystatin	NYS	100		