

Skan

Microplate photometer for fast and standardized reading of all MICRONAUT plates

Technical data

- ▼ Filter: filter wheel with up to 7 filters
(405, 414, 450, 492, 540, 620, 690 nm)
- ▼ Measurement modes: mono-chromatic, bio-chromatic
two point kinetic, multi-chromatic
- ▼ Interfaces: serial and parallel
- ▼ Max. linearity deviation: $\pm 1.5 \%$ at 0-2 A
- ▼ Dynamic range: 0-3.5 A absorbance units
- ▼ Accuracy: - 2 % or 0.007A; typically - 1 %
- ▼ Measures: 140 x 420 x 320 cm

